

# Test Print for your GM3D Printer

Article Number: 13 | Rating: Unrated | Last Updated: Thu, Mar 15, 2018 at 3:14 PM

To verify your printer is calibrated and printing correctly, download the appropriate Test Block file for your specific 3d printer. These were all created to print with **PLA**.

GM3D200 with a 0.4mm Nozzle - GM3D200-4mm-TestBlock.gcode

GM3D200 with a 0.8mm Nozzle - GM3D200-8mm-TestBlock.gcode

GM3D500 with a 0.8mm Nozzle - GM3D500-8mm-TestBlock.gcode

GM3D500/505 R2 with a 0.8mm Nozzle - GM3D500-505R2-8mm-TestBlock.gcode

Posted by: [Shelley Warner](#) - Thu, Mar 15, 2018 at 3:12 PM. This article has been viewed 1346 times.

Online URL: <http://gorillamaker.knowledgebase.co/article.php?id=13>